

Introducing the
Mageleka MagnoMeter™ XRS

A New Instrument for
Particle Characterization and Surface Analysis
of Suspensions and Slurries

What is a *MagnoMeter*?



- ⊗ A device that uses NMR relaxation for non-invasive routine analysis of complex solid-liquid and liquid-liquid formulations
 - ⊗ Samples in water, solvents and melts
 - ⊗ No dilution required, minimal sample preparation
- ⊗ Ideal for applications where speed of measurement and reliability are drivers of performance
 - ⊗ Research & development
 - ⊗ Quality control/quality assurance
 - ⊗ Process control labs

What is NMR relaxation?



- ⊗ NMR relaxation time is a fundamental, intrinsic property of all solids and liquids
 - ⊗ Analogous to intensity of scattered light (particle sizing) and electrophoretic mobility (zeta potential)
- ⊗ NMR relaxation measures how protons react, through their molecular motion, in a magnetic field
- ⊗ Relaxation time of suspensions/slurries is intermediate between that for solid and liquid
 - ⊗ Value depends on specific particle-liquid combination
 - ⊗ Determined directly using an NMR spectrometer

Why use NMR relaxation?



“ In every industrial application, a knowledge and understanding of the molecular structure and dynamics at the particle-liquid interface is critical to improving, or optimizing, suspension and emulsion product performance. ”

NMR Instruments



High resolution NMR

well-known technique for studying molecular structure and identification of compounds



- ⊗ high frequency needed
- ⊗ expensive, complex, sophisticated operation
 - ⊗ intensive training
- ⊗ university and analytical laboratories

Low field NMR

new technique for suspension and emulsion analysis



- ⊗ low frequency optimal
- ⊗ inexpensive, simple benchtop device
 - ⊗ easy operation
- ⊗ industrial R&D, QC/QA and process laboratories

What can you measure?



⊗ Dispersed (wetted) Surface Area of suspensions

⊗ Dispersed particle volume fraction

⊗ Molecular weight (polymers in solution)*

⊗ Polymer and solvent viscosity

⊗ Particle Size (<30nm)

⊗ Relaxation Number (liquids/suspensions/emulsions)

⊗ Kinetic processes

⊗ Adsorption/desorption

⊗ Competitive adsorption

⊗ Colloidal stability

⊗ Presence of para- and ferro-magnetic impurities

⊗ Hydroxyl (OH) number of metal oxides

⊗ Oxygen and water content of solvents

* By calibration

“ The MagnoMeter can be used in an almost unlimited range of applications, and can measure samples at virtually any industrially relevant solids concentration.”

- ⊗ Graphene/Graphene Oxide
- ⊗ Cellulose nanocrystals
- ⊗ Ceramics, refractories
- ⊗ Minerals, metal oxides
- ⊗ Paints and inks
- ⊗ Dyes
- ⊗ Pharmaceuticals
- ⊗ Nanomedicine
- ⊗ Cosmetics
- ⊗ Food emulsions
- ⊗ Agrochemicals
- ⊗ Catalysts
- ⊗ Paper pulp
- ⊗ MOFs

What does the *MagnoMeter* provide?

- ⊗ Direct information about the extent and nature of *any* particle-liquid interface → suspensions and emulsions
 - ⊗ any type of particle, and any liquid – including mixtures
 - ⊗ exceptionally wide concentration (0.01% to 90+%)
 - ⊗ no dilution required, minimal sample preparation
 - ⊗ small samples (typically 0.1 mL; as little as 200 μL)
- ⊗ Complementary information to traditional particle characterization techniques
- ⊗ Intelligence that is not possible with traditional techniques

Advantages of the *MagnoMeter*

General

- ⊗ Any colloidal-size suspension – virtually all liquids
- ⊗ Measure locally with data manipulation at remote terminal
 - ⊗ Ideal for controlled or hazardous environments
- ⊗ Multiple pod sensors facilitate sampling at different locations
- ⊗ Frequency lock
- ⊗ Measure mixture homogeneity

Options and Extensions

- ⊗ Large diameter NMR tube
 - ⊗ viscous liquids, concentrated slurries/emulsions
- ⊗ MRI methods (imaging) in specific samples
- ⊗ Can be adapted for use with an auto-sampler
 - ⊗ multiple sample analysis
- ⊗ Can be adapted for flow-through operation

The *MagnoMeter* finds use from Fundamentals to End-Use Performance

⊗ *Research & Development*

- ⊗ Basic formulation of products: reproducibility, stability (coagulation/flocculation), settling and sedimentation, shelf life

⊗ *Quality Control*

- ⊗ Release of raw materials, batch-to-batch reproducibility of final product

⊗ *Quality Assurance*

- ⊗ Enable release of complex systems which can currently only be characterized by the raw materials used rather than the manufactured product

⊗ *Process Management*

- ⊗ Follow and monitor milling and manufacture processes in almost real-time

Thank you!

For more information, to send samples,
to arrange a demonstration at your
facility, or to speak to a technical
applications specialist, please contact:

info@mageleka.com

+1 617 331 1130

www.mageleka.com